### HI 96821

# Digital Refractometer for Sodium Chloride Measurement for the Food Industry

• Ideal for the analysis of:

Salad Dressings, Cheeses, Condiments, Pickles, Canned Foods, Jarred Foods, Milk, Juices, Energy Drinks, Soups, Brines and Whey

Dual-level LCD

The dual-level LCD displays measurement and temperature readings simultaneously

Automatic Temperature Compensation
For accurate measurements

#### Easy measurement

Place a few drops of the sample in the well and press the READ key

BEPS

(Battery Error Prevention System) alerts the user in the event that low battery power could adversely affect readings.

## IP65 water protection

Built to perform under harsh laboratory and field conditions.

- Quick, accurate results Readings are displayed in approximately 1.5 seconds.
- Single point calibration Calibrate with distilled or deionized water
- Small sample size Sample size can be as small as 2 metric drops.
- Automatic shut-off After three minutes of non-use
- Stainless steel sample well Easy to clean and corrosion resistant
- ABS thermoplastic casing



HANNA offers the HI 96821 digital sodium chloride refractometer to meet the requirements of the food industry. This optical instrument employs the measurement of the refractive index to determine sodium chloride concentration in aqueous solutions used in food preparation. It is not intended for sea water salinity measurements.

The measurement of refractive index is simple and quick and provides the user an accepted method for NaCl analysis. Samples are measured after a simple user calibration with deionized or distilled water. Within seconds the instrument measures the refractive index of the solution. The digital refractometer eliminates the uncertainty associated with mechanical refractometers and is portable for measurements where you need them.

The instrument utilizes internationally recognized references for unit conversion and temperature compensation. It can display the measurement of NaCl concentration 4 different ways: g/100 g, g/100 mL, specific gravity, and °Baumé.

Temperature (in °C or °F) is displayed simultaneously with the measurement (on 3 of the ranges) on the large dual level display along with icons for Low Power and other helpful message codes.



# With Great Products, Come Great Results™

REFRACTOMETERS

# **Easy Operation**

#### Start-up screens

When the HI 96821 is turned on, test screens then the percentage of battery life remaining is shown on the LCD followed by the ready status.

#### Unit selection

SPECIFICATIONS

Range

Resolution

Accuracy

(@20°C/68°F)

g/100 g

g/100 mL

Specific Gravity °Baumé

Temperature

g/100 g

g/100 mL

°Baumé

g/100 g

°Baumé

**Temperature Compensation** 

Minimum Sample Volume

Battery Type / Battery Life

**Dimensions / Weight** 

Measurement Time

Light Source

Sample Cell

Auto-off Enclosure Rating

g/100 mL

Specific Gravity

Temperature

Specific Gravity

Temperature

Just press the RANGE key to cycle through the HI 96821's units of measurement. g/100 g, g/100 mL, Specific Gravity and °Baumé.

Temperature selection can also be easily changed.



#### Calibration

Perform a quick and easy calibration after start-up:

- **1.** Using a plastic pipette, completely cover the prism in the sample well with distilled or deionized water.
- 2. Press the ZERO key

HI 96821

0 to 28

0 to 34 1.000 to 1.216

0 to 26

0 to 80°C (32 to 176°F)

0.1

0.1

0.001

0.1

0.1°C (0.1°F)

±0.2

±0.2

±0.002

±0.2 ±0.3°C (0.5°F)

automatic between 10 and 40°C (50 to 104°F)

approximately 1.5 seconds

100 µL (to cover prism totally)

yellow LED

stainless steel ring and flint glass prism after three minutes of non-use

IP65

9V / approximately 5000 readings

192 x 104 x 69 mm (7.6 x 4.1 x 2.7") / 420 g (14.8 oz.)

#### Measurement

Achieve fast, professional results:

- **1.** Using a plastic pipette, drip sample onto the prism surface until the well is full.
- **2.** Press the READ key and the results are display in the selected units.

## Making a standard sodium chloride solution

To make a Standard NaCl Solution (g/100 g), follow the procedure below:

- Place a container (such as a glass vial or dropper bottle that has a cover) on an analytical balance.
- Tare the balance.
- To make an X NaCl solution weigh out X grams of high purity dried Sodium Chloride (CAS #: 7647-14-5: MW 58.44) directly into the container.
- Add distilled or deionized water to the container so the total weight of the solution is 100 g.

#### Example with g/100 g NaCl:

g/100 g NaCl	10
g NaCl	10.000
g Water	90.000
g Total	100.000



## ORDERING INFORMATION

**HI 96821** is supplied with battery and instruction manual.



# 13

REFRACTOMETERS

www.hannainst.com